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PLEASE READ THESE INSTRUCTIONS CAREFULLY AND KEEP THEM

Low/high elastic orthosis with a thermoformable pad to stabilize the lumbar rachis and relieve pressure

DECLARATION OF CONFORMITY
 As the manufacturer, ORTHOSERVICE AG declares, under its sole responsibility, that this is a class I medical device that has been manufactured and assessed, according to the EU Regulation 2017/745 (MDR). These instructions have been drawn up in application of the aforementioned Regulation. They are intended to ensure an adequate and safe use of the medical device.

TRADEMARKS OF MATERIALS
 Velcro® is a registered trademark of Velcro Industries B.V.

SAFETY PRECAUTIONS
 We recommend that the pressure exerted by the device does not act on parts of the body with wounds, swelling, or welts. It is advisable to not over tighten the device so as to avoid generating areas of excessive local pressure or the compression of underlying nerves and/or blood vessels. Furthermore, we recommend wearing a garment to avoid direct contact with the skin.
 If in doubt on how to apply the device, contact a doctor, physiotherapist, or an orthopedic technician. Read the product composition on the internal label carefully.
 We recommend to not wear the device in the proximity of free flames or strong electromagnetic fields. Do not apply in direct contact with open wounds.

WARNINGS
 It is advisable that the device, designed for the specific indications as set out below, be prescribed by a doctor or a physiotherapist and applied by an orthopedic technician, according to the specific needs of the patient.
 To ensure effectiveness, tolerability, and proper functioning, the application must be carried out with the utmost care. Never alter the adjustment made by the doctor/physiotherapist/ or orthopedic technician. The responsibilities of the manufacturer shall lapse in case of inappropriate use or adaptation. The orthosis is made and to be used by only one patient. Should the orthosis be inappropriately used, the manufacturer declines all responsibility, as provided by the regulation for medical devices. In hypersensitive individuals, direct contact with the skin may cause redness and irritation. In case of pain, swelling, welts, or any other adverse reaction, contact your doctor immediately, and in case of a particularly severe adverse event, report the incident to the manufacturer and the competent authority of your country. The effectiveness of the orthopedic product is only guaranteed when all the components are properly in use.

Item	REF.T1170					
Size	XS	S	M	L	XL	XXL
Circ. waist cm	50/60	60/75	75/90	90/105	105/120	120/135
Front height cm	20	20	20	20	20	20
Back height cm	26	26	26	26	26	26

Colore: black

Item	REF.T1172					
Size	XS	S	M	L	XL	XXL
Circ. waist cm	50/60	60/75	75/90	90/105	105/120	120/135
Front height cm	22	22	22	22	22	22
Back height cm	32	32	32	32	32	32

Colore: black

Option:
 REF. 1078/P Lumbar pad
 REF. 9ST012 Additional reinforced stay 22 cm (Lumbofit70 only) or 26 cm (Lumbofit72 only)

MAINTENANCE
 ☒ Do not bleach ☒ No chemical cleaning
 ☒ Do not iron ☒ Do not tumble-dry
 ☒ Washing instructions: hand wash at max. 30° C with neutral soap. Dry away from heat sources.
 Do not dispose of the product or any of its components into the environment.

INDICATIONS
 • Lumbago
 • Lumbar sciatica and cruralgia
 • Spondylarthrosis and lumbar discopathy
 • Slight traumas to the lumbosacral column and paravertebral lumbar contractures
 • Postoperative recovery
 • **LumboFIT72 only:** Results of vertebral compressions of the lumbar spine
 • **LumboFIT72 only:** Treatment of transverse process fractures

CONTRAINDICATIONS
 Currently no known

CHARACTERISTICS AND MATERIALS
 • Separate, large-scale tension straps provide support and can be individually adjusted
 • Front and side material made of viscose staple fibre and elastomer
 • Back part made of breathable AirX™ fabric: skin-friendly, non-cytotoxic, latex free
 • 4 back linear pockets allow insertion of the optional rigid reinforced stays as a replacement for the ones provided
 • The abdominal section is reinforced with oblique spiral stays
 • Front Velcro® fastening system
 • Adapts to the pathology of the patient

PUTTING ON THE APPLIANCE
 1 Open the orthosis and temporarily fix the Velcro® straps at the lateral straps
 2 Open the orthosis and position with the orange label facing upwards.
 3 Fit the orthosis and fix using the middle Velcro® closure (fig. A).
 4 Tighten the orthosis by fixing the Velcro® straps on the front Velcro® strap (fig. B). ONLY high model (LumboFIT72, REF. T1172): First pull tight the lower and then the upper Velcro® straps (fig. C).
 5 Ensure that the orthosis is exactly in the centre of the body, both on the front as well as on the back (fig. D).

NOTE:
 To ensure optimal stabilisation, mould the rigid stays provided and slide them into the pockets provided in the back.
Attaching the rigid pad:
 - if necessary, heat and mould the plastic part to fit the shape of the patient's body
 - place the padding on the plastic part (concave side-without any writing-fig.L)
 - apply the adhesive Velcro® fasteners centrally on the opposite side to the padding (fig.M)
 - attach the pad to the central Velcro® to the inside of the back of the orthosis

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