LEVITAS®

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ENGLISH

MA LAD662-LAD663 02 A_ENG_09-2020

**DOMUS ALTERNATING PRESSURE
REDISTRIBUTION SYSTEM**

INSTRUCTION MANUAL

INDEX

1. CODES.....	PAG.3
2. INTRODUCTION.....	PAG.3
3. INTENDED USE.....	PAG.3
4. DECLARATION OF CONFORMITY	PAG.3
5. STANDARDS	PAG.4
6. GENERAL WARNINGS	PAG.4
6.1 Electromagnetic compatibility warnings.....	pag. 4
7. SYMBOLS.....	PAG.5
8. KIT CONTENTS.....	PAG.6
9 GENERAL DESCRIPTION	PAG.6
9.1 Pump unit	pag. 6
9.2 Front panel	pag. 6
9.2.1 Pressure Adjust Knob.....	pag. 6
9.2.2 Main Power Switch.....	pag. 6
10. INSTALLATION.....	PAG.7
11. BEFORE USE	PAG.7
12. USING WARNINGS.....	PAG.7
13. OPERATION	PAG.8
13.1 General operation	pag. 8
14. CLEANING	PAG.8
15. STORAGE	PAG.9
16. MAINTENANCE	PAG.9
17. CONDITIONS OF DISPOSAL	PAG.9
18. SPARE PARTS/ACCESSORIES.....	PAG.9
19. TROUBLESHOOTING	PAG.9
20. GUIDANCE AND MANUFACTURER'S DECLARATION	PAG.10
20.1 Guidance and Manufacturer's Declaration - Electromagnetic Emissions.....	pag. 10
20.2 Guidance and Manufacturer's Declaration - Electromagnetic Immunity	pag. 10
20.3 Recommended separation distances between portable and mobile RF communications equipment and this device	pag. 12
21. TECHNICAL FEATURES.....	PAG.13
21.1 Pump	pag. 13
21.2 Mattress.....	pag. 14
22. WARRANTY	PAG.14
23. REPAIR.....	PAG.14
23.1 Warranty repairs	pag. 14
23.2 Repair not covered by warranty.....	pag. 14
23.3 Non-defective device	pag. 15
24. SPARE PARTS	PAG.15
25. EXEMPT CLAUSES	PAG.15

CE I Class Medical Device

REGULATION (EU) 2017/745 OF THE EUROPEAN PARLIAMENT AND OF THE COUNCIL of 5 April 2017 concerning medical devices

1. CODES

- LAD662** Kit mattress and pump with pressure adjustment Domus 1 (bubble pad)
LAD663 Kit mattress and pump with pressure adjustment Domus 2 (single cells pad)

2. INTRODUCTION

Thank you for purchasing a LEVITAS mattress and pump Domus Kit (single cell pad) by Moretti. The pump and mattress kits by Moretti are designed and manufactured to meet all your needs for a practical, correct and safe use. This manual contains small suggestions for proper use of the device you have chosen and valuable advice for your security. You may read the complete manual before using the kit pump and mattress, in case of doubts please contact your dealer, who will help and advise you properly.

NOTE, CAUTION AND WARNING STATEMENTS:

NOTE: Indicate some tips.

CAUTION: Indicate correct operating or maintenance procedures in order to prevent damage to or destruction of the equipment or other property.

WARNING: Calls attention to a potential danger that requires correct procedures or practices in order to prevent personal injury.

3. INTENDED USE

The kit pump and mattress have been designed for the bed sore treatment and prevention, which may occur during long hospital and/or house confinement.

CAUTION!

- Do not use the product for a purpose not indicated in this manual
- Only qualified personnel or trained for the treatment and prevention of bed sore, can operate with the device
- Moretti S.p.A declines all responsibilities for any consequences resulting from an incorrect use of this product and from unauthorized alteration to the frame of the product.
- The manufacturer reserves the right to change the information contained in this document without previous notice.

4. DECLARATION OF CONFORMITY

Moretti S.p.A. declares under its sole responsibility that the product made and traded by Moretti S.p.A. and belonging to the group of DOMUS ANTIDECUPITOS KIT complies with the provisions of the regulation 2017/745 on MEDICAL DEVICES of 5 April 2017.

For this purpose, Moretti S.p.A. guarantees and declares under its sole responsibility what follows:

1. The devices satisfy the requirements of general safety and performance requested by the Annex I of regulation 2017/745 as laid down by the Annex IV of the above mentioned regulation.
2. The devices ARE NOT MEASURING INSTRUMENTS.
3. The devices ARE NOT MADE FOR CLINICAL TESTS.
4. The devices are packed in NON-STERILE BOX.
5. The devices belong to class I in accordance with the provisions of Annex VIII of the above mentioned regulation
6. Moretti S.p.A. provides to the Competent Authorities the technical documentation to prove the conformity to the 2017/745 regulation, for at least 10 years from the last lot production.

Note: Complete product codes, the manufacturer registration code (SRN), the UDI-DI code and any references to used regulations are included in the EU declaration of conformity that Moretti S.p.A. releases and makes available through its channels.

5. STANDARDS

The system has been tested and successfully approved to the following standards:

- EN 60601-1
- EN 60601-1-2
- EN 55011 Class B
- EN 61000-3-2 Class A
- EN 61000-3-3

6. GENERAL WARNINGS

- For proper use of the device, refer to the following manual carefully
- Advise with physician or therapist to determinate the correct adjustment and the correct use of the device
- Keep the packed kit away from heat sources
- SERVICE LIFE- the device use limit is defined by the wear of the parts
- DO NOT allow children to play on or operate the pump
- The user and/or the patient will have to report any serious accident that have occurred related the device to the manufacturer and appropriate authority of the State which the user and/or patient belongs to.

6.1 Electromagnetic compatibility warnings

This equipment has been tested and it complies with medical devices limits belonging to EN 60601-1-2. These limits have been designed to provide reasonable protection against harmful interference in a typical medical installation. This equipment generates, uses and can radiate radio frequency energy and, if not installed and used in accordance with the instructions, may cause harmful interference to other devices in the vicinity. However, there is no guaranteed that the interference will not occur in a particular installation. If this equipment does cause harmful interference to other devices, which can be determined by turning the equipment off and on, the user is encouraged to try to correct the interference following one or more of the suggestions below:

- Reorient or relocate the receiving device
- Increase the separation between the equipment
- Connect the equipment into an outlet on a circuit different from that to which the other devices are connected
- Consult the manufacturer or field service technician for help

7. SYMBOLS

	Product code
	Unique Device Identification
	CE mark
	Manufacturer
	Batch Lot
	Read the instruction manual
	Medical Device
	Conditions of disposal
	Attention, should read the instructions
	Production dates
	Waste Electrical & Electronic Equipment (WEEE) UE/19/2012
	Temperature Limitation
	BF type
	Class II Equipment
IP21 	Protected against solid foreign objects of 12,5 mm and greater. Protection against vertically falling water drops.
	On/Off Button
ON OFF	

8. KIT CONTENTS


The DOMUS series is composed by 2 kind of kits which can satisfy all the requirement about treatment and prevention of complications caused by long hospital and/or house confinement, in particular bedsores.

LAD662 Kit mattress and pump with pressure adjustment Domus 1 (bubble pad)


LAD663 Kit mattress and pump with pressure adjustment Domus 2 (single cells pad)

9 GENERAL DESCRIPTION

9.1 Pump unit

- Front
1. Power Switch
 2. Front Panel

- Rear
3. Air Hose Port
 4. Hooks
 5. Power Cord

9.2 Front panel

9.2.1 Pressure Adjust Knob

Pressure adjust knob controls the air pressure output. Please consult the physician for a suitable setting

9.2.2 Main Power Switch

To turn the pump unit on/off

10. INSTALLATION

Unpack the box and check the package contents for completeness. If there are any damages, please immediately contact your retailer.

1. Place the mattress on a common mattress, with the air hose connectors oriented to the footboard

ATTENTION!

We recommend to place the mattress with interchangeable elements over a mattress with a height of at least 5 cm

2. Hang the pump onto the bed rail (footboard side), or place it on a flat surface
3. Connect air hose connectors from air mattress to the pump unit

ATTENTION!

Ensure the air hoses are not kinked or tucked under mattress

4. Plug the power cord into the electrical outlet

-  NOTE: 1. Make sure the pump unit is suitable for the local power voltage
2. The plug is also served to disconnect the device. Do not position the equipment so that it is difficult to disconnect the device

ATTENTION!

The pump should only be used with mattress recommended by the manufacturer. Do not use it for any other purpose

5. Turn the switch ON

Several installation tips are listed below:

After installation, the extra length of the power cord, if any, should be neatly arranged to avoid any tripping accidents. The EQUIPMENT should be firmly placed at position where users/doctors can access easily.

11. BEFORE USE

Check the correct assembly, paying particular attention to air hoses connector.

12. USING WARNINGS

DANGER - To reduce risk of electrocutione

1. Always unplug this product immediately after using
2. Do not use while bathing
3. Do not place or store this product where it can fall or be pulled into a tub or sink
4. Do not place in or drop into water or other liquid
5. Do not reach for a product that has fallen into water. Unplug immediately

WARNING - To reduce risk of burns, electrocution, fire or injury to people

1. Never leave this product unsupervised when connected to the mains
2. Evaluate patients for entrapment risk according to facility protocol and monitor patient appropriately
3. This system is not intended to be used with patients who have spinal cord injury.
4. Close supervision is necessary when this product is used on or near children.
Electrical burns or shocking accident may occur if a child swallow a small part detached from the device
5. Use this product only for its intended use as described in this manual. Do not use attachment not recommended by the manufacturer
6. Never operate this product if it has a damaged cord or plug, if it is not working properly, or if it has been dropped, damaged, or dropped into water. Return the product to a service centre for examination and repair
7. Keep the cord away from heated surfaces
8. Never block the air openings of this product or place it on a soft surface, such as a bed or couch, where their openings may be blocked. Keep the air openings free of lint, hairs, and other similar particles
9. Never drop or insert any object into any opening or hose
10. Do not modify this equipment without authorization of the manufacturer
11. Mattress cover has passed skin sensitization and skin irritation test. However, if you suspect that you may have had or are having an allergic reaction, please consult a physician immediately
12. Do not leave excessively long tubes in the upper area of the bed, to avoid the possibility of strangulation

13. OPERATION

 NOTE: Always read the operating instruction before use

13.1 General operation

1. Switch on the main power switch on right side of the pump
2. Adjust the pressure setting based upon patient comfort level by turning the pressure adjustment knob clockwise to increase firmness

 NOTE: every time the mattress is set up for use, it is recommended that the pressure first to be set to the max. The user/career can then adjust air mattress weight levels to the desired softness after set up has been completed

14. CLEANING

It is important to follow the cleaning procedures before using the equipment on human bodies.

Wipe the pump unit with a damp cloth and a mild detergent. If other detergent is used, choose one that will have no chemical effects on the surface of the plastic case of the pump unit.

Wipe down the mattress with warm water containing a mild detergent. The cover may also be cleaned by using sodium hypochlorite diluted in water. All parts should be air dried thoroughly before use.

 CAUTION: Do not immerse or soak pump unit.

 CAUTION: Do not use phenolic based product for cleaning.

 CAUTION: After cleaning, dry the mattress without direct exposure sunlight.

WARNING!
Unplug the pump before starting the cleaning procedure
15. STORAGE

For bubble pads:

- Fold the mattress accurately to prevent creases
- Roll from the head end towards the foot end with the valve open

For single cells mattress:

1. Lay the bubble pad or mattress on a flat surface and upside down
2. Roll-up the mattress from the head end towards the foot end
3. Foot-end strap can then be stretched around the rolled pad/ mattress to prevent unrolling

 **CAUTION:** Do not fold, crease or stack the mattresses.

16. MAINTENANCE

The LEVITAS devices by MORETTI are checked carefully and supplied with EC mark, once launched on the market.

The products are intended to offer safe and reliable operation when used or installed according to the instructions provided by Apex Medical. Apex Medical recommends that the system be inspected and serviced by authorized technicians if there are any signs of wear or concerns with device function. Otherwise, service and inspection of the devices generally should not be required.

1. Check main power cord and do not plug it if there is an abrasion or excessive wear.
2. Check mattress cover for signs of wear or damage.
3. Disconnect the air tube from the mattress. Then check the airflow coming from the two air outlets on the pump. They should be alternately delivering air when the pump is set at "alternate" mode.
4. Check the air hoses if there is kink or breaks. For replacement, please local distributor.

17. CONDITIONS OF DISPOSAL

TREATMENT OF BATTERIES (Directive 2012/19/UE).

At the end of its life, the product must not be disposed of along with other domestic waste. The users must dispose of this equipment by bringing it to a specific recycling point for electric and electronic equipment or at retailers that provide this service. By ensuring these batteries are disposed correctly, you will help prevent potentially negative consequences for the environment and

human health which could otherwise be caused by inappropriate waste handling of the battery. To remark the need to dispose of electrical equipment separately, the products is marked with crossed mobile waste bin

18. SPARE PARTS/ACCESSORIES

For spare parts and accessories refer to main catalogue

19. TROUBLESHOOTING

If your questions can't be answered with above information, please contact your local agent directly. They might require a technician to take care of the problem

PROBLEM	SOLUTION
Power is not ON	<ul style="list-style-type: none"> • Check if the plug is connected to mains
The mattress doesn't inflate	<ul style="list-style-type: none"> • Check if the connection between mattress and pump unit is tightly secured • Check if there is any leakage on air cells
Mattress form is loose	<ul style="list-style-type: none"> • Check if all the snap buttons or straps of mattress are all securely fastened • Check if the mattress is fixed to the bed frame by straps
No air produced from some air outlets of the air tube connector	<ul style="list-style-type: none"> • This is normal since there is alternating mode. Air outlets take turns to produce air during their cycle time

⚠ NOTE: If the pressure level is consistently low, check for any leakage (tubes or air hoses). If necessary, replace any damaged tubes or hoses or contact your local qualified dealer for repair.

20. GUIDANCE AND MANUFACTURER'S DECLARATION

20.1 Guidance and Manufacturer's Declaration - Electromagnetic Emissions

This device is intended for use in the electromagnetic environment specified below. The user of this device should make sure it is used in such an environment.

EMISSIONS TEST	COMPLIANCE	ELECTROMAGNETIC ENVIRONMENT - GUIDANCE
RF emissions CISPR 11	Group 1	The device uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment
RF emissions CISPR 11	Class B	The device is suitable for use in all establishments, including domestic establishments and those directly connected to the public low-voltage power supply network
Harmonic emissions IEC61000-3-2	Class A	
Voltage fluctuations / Flicker emissions IEC61000-3-3	Complies	

⚠ Warning:

1. The device should not be used adjacent to or stacked with other equipment. If adjacent or stacked use is necessary, the device should be observed to verify normal operation in the configuration in which it will be used.
2. Use of accessories, transducers and cables other than those specified or provided by the manufacturer of this equipment could result in increased electromagnetic emissions or decreased electromagnetic immunity of this equipment and result in improper operation.
3. Portable RF communications equipment (including peripherals such as antenna cables and external antennas) should be used no closer than 30 cm (12 inches) to any part of the pump, including cables specified by the manufacturer. Otherwise, degradation of the performance of this equipment could result.

20.2 Guidance and Manufacturer's Declaration - Electromagnetic Immunity

This device is intended for use in the electromagnetic environment specified below. The user of this device should make sure it is used in such an environment.

Basic EMC standard	Immunity Test Levels		Compliance Levels	Electromagnetic Environment-Guidance
	Professional healthcare facility environment	HOME HEALTHCARE ENVIRONMENT		
Electrostatic Discharge (ESD) IEC61000-4-2	+8kV contact +15kV air		+8kV contact +15kV air	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30 %.
Electrical fast transient/ burst IEC61000-4-4	+2kV for power supply line +1kV for input/output line		+2kV for power supply line +1kV for input/output line	Mains power quality should be that of atypical commercial or hospital environment
Surge IEC61000-4-5	+ 1 kV line(s) to line(s) + 2 kV line(s) to earth	+ 1 kV line(s) to line(s)	+ 1 kV line(s) to line(s)	Mains power quality should be that of atypical commercial or hospital environment.
Voltage dips, short interruptions and voltage variations on power supply input lines IEC61000-4-11	Voltage Dips: i) 100% reduction for 0.5 period, ii) 100% reduction for 1 period, iii) 30% reduction for 25/30 period, Voltage Interruptions: 100% reduction for 250/300 period		230V	Mains power quality should be that of a typical commercial or hospital environment. If the user of this device requires continued operation during power mains interruptions, it is recommended that the device be powered from an uninterruptible power supply or a battery.
Power frequency (50/60Hz) magnetic field IEC61000-4-8	30 A/m	30 A/m	30 A/m	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.

Basic EMC standard	Immunity Test Levels		Compliance Levels	Electromagnetic Environment-Guidance
	Professional healthcare facility environment	HOME HEALTHCARE ENVIRONMENT		
Conducted RF IEC 61000-4-6	3 Vrms 0,15 MHz - 80 MHz 6 Vrms in ISM bands between 0,15 MHz and 80 MHz 80 % AM at 1 kHz	3 Vrms 0,15 MHz - 80 MHz 6 Vrms in ISM and amateur radio bands between 0,15 MHz and 80 MHz 80 % AM at 1 kHz	6Vrms	Portable and mobile RF communications equipment should be used no closer to any part of this device, including cables, than there commended separation distance calculated from the equation applicable to the frequency of the transmitter. Recommended separation distance $=\sqrt{P} \cdot 150\text{kHz to } 80\text{MHz}$ $=0.6\sqrt{P} \cdot 80\text{MHz to } 800\text{MHz}$ $=1.2\sqrt{P} \cdot 800\text{MHz to } 2.7\text{GHz}$ Where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and d is the recommended separation distance in meters (m). Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey, should be less than the compliance level in each frequency ranged. Interference may occur in the vicinity of equipment marked with the following symbol: 
Radiated RF EM Fields IEC61000-4-3	3 V/m 80 MHz to 2.7 GHz 80 % AM at 1 kHz 385-6000 MHz, 9-28V/m, 80% AM(1kHz) pulse mode and other modulation	10 V/m 80 MHz to 2.7 GHz 80 % AM at 1 kHz 385-6000 MHz, 9-28V/m, 80% AM(1kHz) pulse mode and other modulation	10V/m	
NOTE 1: UT is the a.c. mains voltage prior to the application of the test level NOTE2: At 80 MHz and 800 MHz, the higher frequency range applies. NOTE 2: These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people				
a)Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the device is used exceeds the applicable RF compliance level above, the device should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as reorienting or relocating the device. b) Over the frequency range 150 kHz to 80 MHz, field strengths should be less than 10 V/m.				

20.3 Recommended separation distances between portable and mobile RF communications equipment and this device

This device is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer or the user of this device can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and this device as recommended below, according to the maximum output power of the communications equipment.

Rated maximum output power of transmitter W	Separation distance according to frequency of transmitter m		
	150 kHz to 80 MHz =√P	80 MHz to 800 MHz =0,6√P	800 MHz to 2,7 GHz d=1,2√P
0.01	0.1	0.06	0.12
0.1	0.31	0.19	0.38
1	1	0.6	1.2
10	3.1	1.9	3.8
Rated maximum output power of transmitter W	Separation distance according to frequency of transmitter m		
100	10	6	12

21. TECHNICAL FEATURES

21.1 Pump

MODEL	LAD662	LAD663
Power supply	230V/50Hz	230V/50Hz
Fuse T1AL	250V	250V
Power cord	2 mt	2 mt
Pressure output	30~80mmHg	30~80mmHg
Air release	>4LPM	>4LPM
Cycle time	9.6 min	9.6 min
Dimensions	25x12,5x8,5 cm	25x12,5x8,5 cm
Weight	1,4 Kg	1,4 Kg
Environment temperature	Operation: +10°C/+40°C Storage: -15°C/+50°C Shipping: -15°C/+70°C	Operation: +10°C/+40°C Storage: -15°C/+50°C Shipping: -15°C/+70°C
Environment humidity	Operation: 10%/90% non condensing Storage: 10%/90% non condensing Shipping: 10%/90% non condensing	Operation: 10%/90% non condensing Storage: 10%/90% non condensing Shipping: 10%/90% non condensing
Atmospheric pressure	700/1013,25 hPa	700/1013,25 hPa
Classification	Class II, Type BF, IP21 Applied part: air mattress Not suitable for use in the presence of a flammable anesthetic mixture (No AP or APG protection)	Class II, Type BF, IP21 Applied part: air mattress Not suitable for use in the presence of a flammable anesthetic mixture (No AP or APG protection)

21.2 Mattress

MATTRESS	MATTRESS BUBBLE PAD (LAD662)	ELEMENT MATTRESS (LAD663)
Model	6,4 cm bubble pad	10,2 cm Overlay mattress
Components	130 bubbles	18 cells
Material	Medical-grade PVC	Cover: Nylon/PU Cell: Nylon/PVC
MATTRESS	MATTRESS BUBBLE PAD (LAD662)	ELEMENT MATTRESS (LAD663)
Dimension	196x90x6,4 cm	190x80x10,2 cm
Weight	2,3 Kg	3,8 Kg
Max user's weight	100 Kg	140 Kg
Thickness	0,3	-

22. WARRANTY

Moretti products are guaranteed from material or manufacturing faults for a period of (two) 2 years from the purchasing date, subject to any limitations and exclusions specified below. The warranty does not apply in case of damages caused by improper use, abuse or modification of the product, and the warranty will not be valid if the instructions for use are not strictly followed. The correct intended use of the device is specified in this manual.

Moretti is not responsible for damages resulting from personal injury or otherwise caused by or in relation to wrong installation or improper use. Moretti does not guarantee the products against damage or defects in the following conditions: natural disaster, not authorized maintenance or repair, damage due to power supply problems (where provided) use of spare parts not supplied by Moretti, improper use, not authorized alteration, shipping damage (other than original shipment by Moretti), or in case of insufficient maintenance as indicated in the manual. The warranty doesn't cover components subject to wear and tear during the correct use of the device.

23. REPAIR

23.1 Warranty repairs

In the event that a product Moretti to be defective in material or manufacturing faults during the warranty period, Moretti will confirm with customer if the fault can be covered from warranty. Moretti, in its sole discretion, can repair or replace the item under warranty at a Moretti dealer or at Moretti location . The labor cost for the repair may be charged to Moretti if the repair is covered by warranty. Repair or replacement does not renew or extend the warranty

23.2 Repair not covered by warranty

A product not covered by warranty can be returned for repairing only after receiving prior authorization from Moretti customer service. The labor and shipping costs for a repair not covered by warranty will be fully borne by the customer or dealer. Repairs of products not covered by warranty are warranted for six (6) months from the date you receive the product repaired.

23.3 Non-defective device

The customer will be informed if, after having examined and tested returned product, Moretti claims that the product is not defective. In this case the good will be returned to the customer, all shipping costs will be paid by customer.

24. SPARE PARTS

For spare parts only refer to MORETTI master catalogue.

25. EXEMPT CLAUSES

Except as expressly specified in this warranty and within the limits of the law, Moretti does not offer any representation, warranty or condition, express or implied, including any representations, warranties or conditions of merchantability, fitness for a particular purpose, non-infringement and interference. Moretti does not warrant that the use of the product Moretti will be uninterrupted or error free. The duration of any implied warranties that may be imposed by law are limited to the warranty period, within the limits of the law. Some states do not allow limitations on how long an implied warranty lasts or the exclusion or limitation of incidental or consequential damages in connection with products for consumers. In such states or countries, certain disclaimers and limitations of this warranty may not apply to you. This warranty is subject to change without notice.

**WARRANTY CERTIFICATE**

Product _____

Purchased on (date) _____

Retailer _____

Address _____ **Town/city** _____

Sold to _____

Address _____ **Town/city** _____

Via Bruxelles, 3 - Meleto 52022 Cavriglia (Arezzo) Tel. +39 055 96 21 11

www.morettispa.com email: info@morettispa.com

MADE IN P.R.C.

MORETTI S.P.A.

Via Bruxelles, 3 - Meleto
52022 Cavriglia (Arezzo)

Tel. +39 055 96 21 11
Fax. +39 055 96 21 200

www.morettispa.com
info@morettispa.com