polfit 17 / 19 / 21

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ENGLISH

Please read the instructions carefully and thoroughly

DECLARATION OF CONFORMITY

As the manufacturer, ORTHOSERVICE AG declares with sole responsibility that this product is a Class I medical product and has been manufactured in accordance with the basic requirements of the directives 93/42/EEC and 2007/47/EEC. These operating instructions were created in accordance with the specification of the above-mentioned directives. They aim to ensure the safe and correct use of the medical product.

TRADEMARKS OF MATERIALS

Velcro® is a registered trademark of Velcro Industries B.V.

SELECTION/SIZE

Size	S	M	L	XL
Wrist circumf. cm	13/16	16/18	18/20	20/23
Length cm (polfit 17 + 19)	20	20	20	20
Length cm (polfit 21)	26	26	26	26

Colour: black - specify right or left

MAINTENANCE

- hand wash at max. 30° C no chemical cleaning
- do not iron do not tumble-dry
- do not bleach one patient only

- The product is maintenance-free.
- Washing instructions: Remove the splint. Hand wash in lukewarm water with mild soap (at max. 30° C); rinse thoroughly. Dry away from heat sources.
- Do not dispose of the device or any of its components.

PRECAUTIONS

It is advised that the pressure exerted by the product does not act on parts of the body suffering from wounds, swelling or lumps. If in doubt about how to apply, contact an orthopedic technician. The device should not be worn in the vicinity of open flames or strong electromagnetic fields.

WARNINGS

The product should be prescribed, used under medical supervision and applied by an orthopedic technician, the competent authority figure for the application as well as for information on the safe use, according to individual needs. To ensure the effectiveness, tolerability, and correct operation, it is necessary that the device is put on with the utmost care. Any changes to the structure or adjustment of the device must be prescribed by a physician and performed by an orthopedic technician. Under no circumstances should the adjustment made by the physician/orthopedic technician be altered. Use of the product is recommended for a single patient. With hypersensitive patients, direct contact with skin may cause redness or irritation. In the case of onset of pain, swelling, lumps or any unusual reaction, immediately contact your physician.

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ENGLISH

CHARACTERISTICS AND MATERIALS
POLFIT 17
 • Aluminum thumb splint preformed and malleable (removable)
 • Aluminum palm splint preformed and malleable (removable)
 • 3D breathable fabric structure
 • Reinforcing dorsal plastic stays
 • Fastening with buckles and Velcro® straps

POLFIT 19
 • Aluminium palm splint, preformed and malleable (removable)
 • 3D breathable fabric structure
 • Reinforcing dorsal plastic stays
 • Fastening with buckles and Velcro® straps

POLFIT 21
 • Full back opening
 • 3D breathable fabric structure
 • Closure with four rings and Velcro® straps
 • Dorsal reinforcing rods in plastic
 • Rigid palmar rod in pre-shaped and mouldable (removable) aluminium

INDICATIONS

POLFIT 17
 • Wrist distortion
 • Wrist tendonitis
 • Arthritis and arthrosis
 • Consequences of wrist fracture or after early removal of the plaster
 • Arthrosis of trapezio-metacarpal joint and metacarpophalangea of the thumb
 • Damage to ulnar section of metacarpophalangea of the thumb

POLFIT 19
 • Wrist distortion
 • Compound fracture of the forearm of radius and ulna
 • Wrist tendonitis
 • Arthritis and arthrosis
 • Consequences of fracture of the wrist or after early removal of the plaster

POLFIT 21
 • Treatment of compound fractures of the distal mega-epiphysis of the radius and ulna that the orthopaedist deems stable
 • Post-surgical treatment of unstable distal fracture of the radius and ulna
 • Wrist distortion
 • Tendonitis
 • Fractures of the carpal bone (except the scaphoid) whether treated surgically or conservatively
 • Complications of rheumatoid arthritis
 • Post-surgical treatment of metacarpal fractures
 • Radio-carpal arthrosis

CONTRAINDICATIONS POLFIT 17 / 19 / 21

No known

PUTTING ON THE APPLIANCE POLFIT 17

ADAPTING FOR THE DOCTOR / TECHNICIAN

1 Adjust the palm metal splint to adapt it to the patient's wrist after removing it from the brace. Put the splint back into the brace (Fig. 1). Repeat the procedure for the thumb metal splint.

PUTTING ON THE APPLIANCE FOR THE PATIENT

- Open all the Velcro® straps and pull the straps out of the buckle (fastening them provisionally on themselves).
- Apply the brace such that the metal splint is positioned at the palm and the thumb is inserted into the provided support (Fig. A). Close the fabric insert by using the Velcro® fastening (Fig. B).
- In case a complete opening of the brace is not necessary, put the hand under the above mentioned precautions (positioning the metal splint at the palm) into the closed brace (Fig. C).
- Fixing the brace:
 - Wrap the central strap around the wrist, pass it through the relevant buckle, tighten it and fasten it with the Velcro® strap (Fig. D).
 - Pass the remaining straps through the relevant buckles, tighten them and fasten them with the Velcro®, first the hand one (Fig. E) and then the forearm one (Fig. F).

5 After fixing the brace fasten the Velcro® strap of the thumb opening (Fig. G).

PUTTING ON THE APPLIANCE POLFIT 19 + POLFIT 21

ADAPTING FOR THE DOCTOR / TECHNICIAN

1 Adjust the palm metal splint to adapt it to the patient's wrist after removing it from the brace. Put the splint back into the brace (Fig. 1).

PUTTING ON THE APPLIANCE FOR THE PATIENT

- Open all the Velcro® straps and pull the straps out of the buckle (fastening them provisionally on themselves).
- Apply the brace such that the metal splint is positioned at the palm and the thumb is inserted into the provided opening (Fig. A). Close the fabric insert by using the Velcro® fastening (Fig. B).
- In case a complete opening of the brace is not necessary, put the hand under the above mentioned precautions (positioning the metal splint at the palm) into the closed brace (Fig. C).
- Fixing the brace:
 - Wrap the central strap around the wrist, pass it through the relevant buckle, tighten it and fasten it with the Velcro® strap (Fig. D).
 - Pass the remaining straps through the relevant buckles, tighten them and fasten them with the Velcro®, first the hand one (Fig. E) and then the forearm one (Fig. F).

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MAINTENANCE

- Do not bleach No chemical cleaning
- Do not iron Do not tumble-dry
- Washing instructions: Hand wash in lukewarm water with mild soap; rinse well. Leave to dry away from direct heat.

Do not dispose of the product or any of its components into the environment.

PRECAUTIONS

It is advised that the pressure exerted by the product does not act on parts of the body suffering from wounds, swelling or lumps. The product should not be tightened too firmly so as to avoid the development of pressure points or compression of the nerves or blood vessels underneath. It is advisable to wear the product over a garment, avoiding direct contact with the skin. If in doubt about how to apply, contact an orthopaedic technician. The device should not be worn in the vicinity of open flames or strong electromagnetic fields.

WARNINGS

The product may only be used in the event of medical indication or when prescribed by a physician and must be applied by an orthopaedic technician, the competent authority figure for the application as well as for information on the safe use, according to individual needs. To ensure the effectiveness, tolerability, and correct operation, it is necessary that the device is put on with the utmost care. Any changes to the structure or adjustment of the device must be prescribed by a physician and performed by an orthopaedic technician. Under no circumstances should the adjustment made by the physician/orthopaedic technician be altered. Use of the product is recommended for a single patient. With hypersensitive patients, direct contact with skin may cause redness or irritation. In the case of onset of pain, swelling, lumps or any unusual reaction, immediately contact your physician.

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ENGLISH

CHARACTERISTICS AND MATERIALS

- Two thumb splints
- Cotton-blend fabric support with padding around the splints
- Adjustable Velcro® thumb strap
- Velcro® wrist strap

INDICATIONS

- Wrist distortion
- Wrist tendonitis
- Arthritis and arthrosis
- Consequences of wrist fracture or after early removal of the plaster
- Arthrosis of trapezio-metacarpal joint and metacarpophalangea of the thumb
- Damage to ulnar section of metacarpophalangea of the thumb

CONTRAINDICATIONS

Currently no known

PUTTING ON THE APPLIANCE

- 1 Slide the thumb into the support (fig. A).
- 2 Adjust the thumb strap to fit and fasten with the Velcro® strap (fig. B).
- 3 Fasten the Velcro® wrist strap (fig. C).
Next place the cuff (ref.PS500 or ref.PS501) over the thumb support.
- 4 Undo the Velcro® straps (fasten them back on themselves for the moment).
- 5 Position the wrist so that the metal stay is on the same side as the palm of the hand and the thumb and thumb support pass through the hole provided (fig. D1/D2).
- 6 Fasten the brace:
 - Fasten the dorsal Velcro® strap (fig. E).
 - Pass the straps through the rings provided, tighten and fasten with Velcro® (fig. F).

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