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## PLEASE READ THESE INSTRUCTIONS CAREFULLY AND SAVE THEM

### Ankle foot orthosis

#### DECLARATION OF CONFORMITY

As the manufacturer, ORTHOSERVICE AG states under its sole responsibility that this product is a Class I medical device that has been manufactured in accordance with the requirements of Directives 93/42/EEC and 2007/47/CE. These instructions have been written in conformity with these directives, and their aim is to ensure the safe and proper use of this medical device.

#### TRADEMARKS OF MATERIALS

Velcro® is a registered trademark of Velcro Industries B.V.

#### PRECAUTIONS

It is advised that the pressure exerted by the product does not act on parts of the body suffering from wounds, swelling or lumps. The product should not be tightened too firmly so as to avoid the development of pressure points or compression of the nerves or blood vessels underneath. It is advisable to wear the product over a garment, avoiding direct contact with the skin. If in doubt about how to apply, contact an orthopaedic technician. Please read the composition of the product on the inside label carefully. The device should not be worn in the vicinity of open flames or strong electromagnetic fields.

#### WARNINGS

The product has been designed for the specific indications stated below and must be prescribed by a doctor and fitted by an orthopaedic technician. The orthopaedic technician is responsible for fitting the product and providing information about how to use it safely in accordance with individual requirements. The product must be fitted with the utmost care to ensure its efficacy, tolerance and correct function. Any changes to the structure or adjustment of the device must be prescribed by a physician and performed by an orthopaedic technician. Under no circumstances should the adjustment made by the physician/orthopaedic technician be altered. The manufacturer is no longer liable if the product is used or adapted inappropriately. The orthosis is intended to be used only by a single patient; if it is used for more than one patient, the manufacturer accepts no responsibility or liability in accordance with medical device legislation. Direct contact with the skin may cause reddening or irritation in hypersensitive individuals. In the case of onset of pain, swelling, lumps or any unusual reaction, immediately contact your physician. The orthopaedic efficacy of the product may only be guaranteed when all of its components are used.

#### WARRANTY

This warranty is governed by the statutory provisions of the country in which the product was purchased. Where relevant information concerning warranties and the relationship between the seller and buyer shall apply. If you believe that the terms of the warranty may be applicable, first contact the retailer where the product was purchased. Orthoservice will repair or replace all or part of the device and its accessories if defects in materials and workmanship arise during the period of six months from the date of purchase.

#### SELECTION/SIZE

Code	Size	REF.20050		
Measures	XS	S	M	L
Ankle circumference cm	13/17	17/21	21/25	25/29

Colour: black - fits right and left

#### SELECTION/SIZE ACCESSORY

Code	Size	REF.20051		
Measures	S	M	L	
Midfoot circumference cm	17/21	21/25	25/29	

Colour: black - specify right or left

#### MAINTENANCE

- Do not bleach
- No chemical cleaning
- Do not iron
- Do not tumble-dry

Washing instructions: hand wash at max. 30 °C with neutral soap.

Dry away from heat sources.

Do not dispose of the product or any of its components into the environment.

#### INDICATIONS

- Flaccid paralysis of the foot
- Light to moderate spastic equinus foot

#### CONTRAINDICATIONS

Currently no known

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- To be worn ONLY with lace-up shoes.
- Low profile.
- Breathable micropore fabric.
- A soft pad attaches the support to the shoe.
- Tensioning elastic regulates the degree of flexion.
- Gel padding over the Achilles tendon.
- Fastens with Velcro® straps.
- Additional midfoot band for wearing with bare feet (optional).

#### PUTTING ON THE APPLIANCE

- A Separate the three parts of the support (shoe pad, tensioning elastic and the tibial cuff) and fully open the tibial cuff.
- B With the tibial cuff open, rest it above the back of the ankle, with the gel padding over the Achilles tendon and the bend upwards (fig.1a), and wrap it around the ankle using the two Velcro® straps (fig.1b). Leave the D-ring completely free and uncovered.
- C Insert the shoe pad between the laces and the tongue of the shoe (fig.2).
- D Put on the shoe and lace it up (fig.3).
- E With the tensioning elastic folded in two, pass it through the D-ring from top to bottom (fig.4).
- F Attach the folded end of the elastic to the shoe pad (fig.5).
- G Finally, attach the elastic to the support and achieve the desired degree of tension by fastening the two Velcro® ends level with the malleoli (fig.6).

#### ACCESSORY: MIDFOOT BAND

##### CHARACTERISTICS AND MATERIALS

- To be worn with bare feet or open shoes together with the Affix.
- Low profile.
- Soft, breathable micropore fabric.
- Anatomically designed to fit the left or right foot.
- Fastens with Velcro®.

#### Putting on the appliance

- A With the tibial cuff open, rest it above the back of the ankle, with the gel padding over the Achilles tendon and the bend upwards (fig.10a), and wrap it around the ankle using the two Velcro® straps (fig.10b). Leave the D-ring completely free and uncovered.
- B Fit the band onto the midfoot with the ring on the top (fig.11).
- C With the tensioning elastic folded in two, pass it through the D-ring from top to bottom (fig.12).
- D Attach the folded end of the elastic to the ring on the midfoot band (fig.13).
- E Finally, attach the elastic to the support and achieve the desired degree of tension by fastening the two Velcro® ends level with the malleoli (fig.14).