Description

CliniSorb is a synthetic absorbable braided sterile surgical suture composed of homopolymers of glycolide (100%). CliniSorb is dyed violet with D & C violet No.2. CliniSorb meets all the requirements established by the B.P., U.S.P. and E.P.

Indications

CliniSorb is indicated for soft tissue approximation and/or ligation, including use in ophthalmic procedures, but not for use in cardiovascular and neurological procedures.

Actions

CliniSorb elicits a minimal acute inflammatory reaction in tissues, which is followed by gradual encapsulation of the suture by fibrous connective tissue. Progressive loss of tensile strength and eventual absorption of CliniSorb Synthetic Absorbable Sutures occurs by means of hydrolysis, where the polymer degrades to glycolic acid which is subsequently absorbed and metabolized by the body. Absorption begins as a loss of tensile strength without appreciable loss of mass. Implantation studies in animals indicate that CliniSorb retains approximately 50% of its tensile strength at 2 weeks post-implantation, with approximately 20% remaining at three weeks. Absorption of Absorbable Synthetic Suture is essentially complete between 60 and 90 days.

Contra-Indications

This suture, being absorbable, should not be used where extended approximation of tissue is required.

Warning

As with any foreign body, prolonged contact of this or any other suture with salt solutions such as those found in the urinary tract, may result in calculus formation. Users should be familiar with surgical procedures and techniques involving absorbable sutures before employing CliniSorb synthetic absorbable suture for wound closure, as a risk of wound dehiscence may vary with the site of application and suture material used. Acceptable surgical practice should be followed with respect to drainage and closure of contaminated or infected wounds. The use of this suture may be inappropriate in patients suffering from conditions that may delay wound healing. As this is an absorbable suture material, the use of supplemental non absorbable sutures should be considered by the surgeon in the closure of sites undergoing expansion, stretching or distension, or which may require additional support.

Precautions

Under some circumstances, notably orthopaedic procedures, immobilisation by external support may be employed at the discretion of the surgeon. Skin sutures which must remain in place longer than seven days may cause localised irritation and should be snipped off or removed as indicated. In handling this or any other suture material care should be taken to avoid damage. Avoid crushing or crimping damage due to the application of surgical instruments, such as forceps or needle holders. Adequate knot security requires the accepted practice of flat, square ties, with additional throws as warranted by the surgical circumstances and experience of the surgeon.

Sterility

CliniSorb sutures are sterilized by Ethylene oxide gas. Do not use if package is opened or damaged. Discard opened unused sutures.

Storage

Recommended storage conditions: between 15°C to 30°C, away from direct heat. Do not use after expiry date.

Adverse Reactions

Adverse reactions associated with the use of this device include wound dehiscence, failure to provide wound support in closure of the abdomen, chest joints, and other sites where expansion, stretching or distension occur, failure to provide adequate wound support in elderly, malnourished or debilitated patients, or in patients suffering from cancer, anaemia, obesity, diabetes, infected or other conditions which may delay wound healing, enhance bacteria infectivity, minimal acute inflammatory tissue reaction, localised irritation when skin closures are left in place for greater than seven days, calculi formation in urinary and biliary tracts when prolonged contact with salt solutions, such as urine and bile occurs, and pain and oedema at the wound site

How supplied

CliniSorb Synthetic Absorbable Braided Suture is available in sizes 6/0 through 2, undyed or dyed violet with D&C violet No.2 and coated. The suture is supplied in non-needled standard lengths (ligatures) or attached to various needle types in boxes of one dozen.

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CLINISUT (PTY) LTD