

Barella corsia

Ward stretcher

K101

Barella da corsia su carrello

Ward stretcher on chassis

Caratteristiche costruttive

Telaio piano rete - realizzato in tubo d'acciaio a sezione ovale mm. 50x25x1,5 carrello inferiore realizzato in tubo d'acciaio a sezione tonda diametro differenziato.

Maniglia di tiro/spinta - in tubo di acciaio sezione tonda Ø 32x1,5 mm.

Finitura - completamente verniciata con polveri epossidiche atossiche a 200°. Colore standard Grigio RAL 7035.

Piano rete - a due sezioni in tondino di acciaio elettrosaldato Ø 5 mm. a maglie rettangolari 100x50 mm. zincato ed estraibile.

Sezione schienale - regolabile a mezzo di una cremagliera manuale a 7 posizioni. Sezione gambe fissa.

Ruote - 4 ruote girevoli in gomma sintetica Ø 200 mm. di cui due dotate di freno indipendente.

Paracolpi - ai quattro angoli in resina antiacida, antiurto.

Barella con carrello - In caso di necessità la parte superiore della barella può essere asportabile per l'inserimento in autolettiga o per altre particolari necessità.

Technical specifications

Platform - All-welded mainframe made of different sections steel tube. Sections frame made of oval steel tube 50x25x1.5 mm.

Push/Pull handle - Made of steel tube Ø 32x1,5 mm.

Surface finish - Epoxy gloss finish standard color Grey RAL 7035.

Mattress platform section - Two section mattress platform fitted with removable zinc-plated grids made of Ø 5 mm. steel rods and rectangular mesh 100x50 mm.

Backrest - Adjustable by means of a manual 7 positions rack device. Fixed leg section.

Castors - Four castors Ø 200 mm. made of synthetic rubber two of which with independent brake.


Bumpers - Four bumpers at the corners made of anti-acid and anti-shock resin.

Ward stretcher with trolley - If necessary the upper part is detachable to insert in ambulance or other use.



KSP®

tutto più semplice



Accessori

K 9101 Asta portaflebo in acciaio cromato a 2 ganci in ABS.

K 9102 Coppia spondine in acciaio verniciato con dispositivo di bloccaggio automatico.

K 9103 Materassino in sky ignifugo.

K 9104 Cestello in acciaio verniciato.

K 9105 Porta bombole in acciaio cromato.

Accessories

K 9101 Infusion pole made of chrome-plated steel tube with two plastic hocks.

K 9102 Couple of epoxy gloss finish side-rails with automatic locking device.

K 9103 Fire-proof sky mattress.

K 9104 Epoxy gloss finish steel lineum holder basket.

K 9105 Chrome plated steel gas bottle-holder.

