




➔ Articolata confortevole  
sedile basculante



Larghezza braccioli variabile  
Adjustable width between  
armrests



**789 Poltrona articolata basculante**

**Standard:** Super imbottiture anatomiche. Tessuti e imbottiture ignifughi. Schienale regolabile e facilmente ribaltabile sul sedile.  
**Sedile basculante manualmente.** Cuscinetti appoggiamambe anatomici integrali, girevoli lateralmente e regolabili in altezza. Fianchi estraibili. Fianchetti paravesti. Freni. Gomme soft. W.C. anteriore igienizzabile. Maniglie singole di spinta. Ruote antiribaltamento.  
**Optional:** Prolunga. Tavolino.

**789 Pivoting, articulated armchair**

**Standard:** Super anatomical paddings. Fire-retardant paddings and upholsteries. Adjustable backrest that folds onto seat. **Manual pivoting seatrest.** Swinging and adjustable in height legrests with anatomical paddings. Removable arms. Side protections for clothes. Brakes. Soft wheels. Easily cleanable W.C. open in front. Push handles. Anti-overturn small wheels.  
**Optional:** Headrest. Small table.

**Classificazione ISO 12.21.06.003**

