






➔ Supera ogni ragionevole ostacolo



-  Kg 12
- +
-  Kg 3,5
- +
-  Kg 1,5
- +
-  Kg 35
- +
-  + Kg 33 BATTERY

**Kg 85** =

**Every Way** Sedia elettronica basculante per esterno smontabile per trasporto in auto

**Standard:** Every Way con i potenti motori di 500 watt e le grandi ruote è pura espressione di potenza, inarrestabile su ogni terreno. Freni elettromagnetici. Batterie 130 A senza manutenzione. Velocità 8 km/h. Autonomia 40 km. Sedile facilmente basculante. Imbottiture e tessuti ignifughi. Imbottiture super anatomiche. Dotata di selezionatore 5 velocità e fanaleria regolamentare. Carica batterie. Cintura di sicurezza. Ruote antiribaltamento. Freni elettromagnetici.














**Optional:** Comando retrattile. Appoggiagambe regolabili in altezza. Sedile e schienale regolabili elettricamente.

**Every Way** Electronic rocking chair for external use disassembled for easy transport on car

**Standard:** Every Way with strong 500 watt motors and big wheels, it's pure expression of power, unstoppable on every soil. Electromagnetic brakes. 130 A batteries without maintenance. Speed: 8 km/h - Range: 40 km. Easily rocking seat. Padding and upholsteries fire-retardant. Super anatomical padding. Equipped with 5 gears selector and regulation lights with indicators. Electromagnetic brakes.

**Optional:** Retractable control box. Adjustable in height legrests. Electrically adjustable seat and backrest.

Classificazione ISO 12.21.27.009+12.24.03.163

												
cm Ø 38	cm Ø 26	L 40 - 45	L 65	L 80 - 105	h 48 - 58	Standard	20%	Portata	150 kg			