

➔ Classica e bella per ogni uso
Larghezza esterna solo 56 cm



Nuovo bracciolo
brevetto Surace
New armrest design
Surace patented



Schienale piegato di serie
Standard folding backrest



+ Kg 28
BATTERY

710 Magic Sedia elettronica pieghevole per interno, esterno, viaggio

Standard: Schienale pieghevole. Motori a presa diretta 320 watt. Regolatore di velocità. Batterie 100 A senza manutenzione. Autonomia 6 h (40 km). Velocità 7 km/h. Carica batterie. Freni elettromagnetici. Fianchi desk girevoli. Pedane girevoli ed estraibili. Cintura di sicurezza. Ruote antiribaltamento. Comando destro o sinistro. Freni parcheggio.

Optional: Comando allo schienale. Appoggiategambe regolabili girevoli ed estraibili.

710 Magic Electronic folding chair for internal, external use and travel

Standard: Folding backrest. Motors: 320 watt with direct tap. Speed governor. Batteries: 100 A without maintenance. Range: 6 h (40 km). Speed: 7 km/h. Battery charger. Electromagnetic brakes. Desk and tipping armrests. Swinging and detachable footrests. Chest strap. Anti-overnturn small wheels. Right or left control box. Parking brakes.

Optional: Electric control for backrest. Elevating, swinging and detachable legrests.

Classificazione ISO 12.21.27.006

cm Ø 33	cm Ø 21	L 44	L 56	L 75 - 98	h 52	Standard	20%	Optional	Portata