




Agile su strade impegnative



Larghezza braccioli
variabile
Adjustable width



+ Kg 54
BATTERY

750 Maxi Sedia elettronica per uso esterno per percorsi impegnativi

Standard: Con i potenti motori da 700 watt, permette di percorrere ogni tipo di terreno e pendenza. Freni elettromagnetici. Batterie 200 A senza manutenzione. Autonomia 6 h (40 km). Velocità 10 km/h. 5 velocità. Carica batterie. Fanaleria regolamentare con indicatori di direzione. Ammortizzatori. Schienale regolabile. Braccioli regolabili in altezza. Super imbottiture ignifughe. Cintura di sicurezza. Comando destro o sinistro.

Optional: Schienale regolabile elettricamente. Cestino schienale.

750 Maxi Electronic chair for external use for difficult tracks

Standard: With its strong 700 watt motors, allows to drive every type of soil or gradient. Electromagnetic brakes. Batteries: 200 A without maintenance. Range 6 h (40 km). Speed: 10 km/h. 5 gears. Battery charger. Regulations lights with indicators. Shock absorbers. Adjustable backrest. Adjustable armrests. Fire-retardant paddings and upholsteries. Chest strap. Right or left control box.

Optional: Electrically adjustable backrest. Back basket.

Classificazione ISO 12.21.27.009+12.24.03.163

cm Ø 38	cm Ø 25	L 43 - 48	L 48 - 53	L 66 - 71	L 84 - 120	h 55	Standard			20%	Portata 150 kg