



Air 35 Surace: solo 35 cm di ingombro massimo. Utile e svelta nei passaggi stretti

Air 35 Surace: only 35 cm of width. Useful and quick for narrow passages



Air35 Sedia per aerei, aeroporti, stazioni ferroviarie, ferrovie, società di navigazione, croce rossa, ambulanza...

Standard: Telaio tubo acciaio cromato molto robusto, prolunga schienale facilmente piegabile, 2 Ruote diam. 210 mm soft, 2 Cinture di sicurezza allo schienale, 1 Cintura di sicurezza al sedile, fascia ai talloni, tessuto ignifugo.

Optional: 1) Freni a tamburo, 2) 2 Rotelline anteriori diam. 80, 3) Cassettino sotto sedile porta documenti, 4) Logo cliente retro schienale

UK Air35 Wheelchair for airplanes, airports, railway stations, railroads, navigation companies, red cross ambulances...

Standard: Chromed steel tube frame very strong, easily foldable headrest, 2 Diam. 210 mm soft wheels, 2 Security belts at the backrest, 1 Security belt at the seat, heels stripe, fireproof fabric.

Optional: 1) Drum brakes, 2) 2 diam. 80 mm front castors, 3) Storage drawer under the seat, 4) Customer's logo at backrest's rear.

