



Respiratory Therapy

XT-III

CPAP

XT series, the new member has just joined APEX family to care for the health of yours and the relief of your OSA symptom. The lightweight, portable, and streamlined design of XT series is for the convenience and comfort of yours on the run touring around the world wherever you want it to be.

XT series is introduced by APEX with the best interest of yours in mind and for someone who is close to you so they can be at ease knowing you are in good hands with proper care. Please, don't leave home without your XT series.





Features:

- Advanced model with superior features
- Enhanced pressure sensor for better pressure accuracy
- Automatic leak adjustment and manual altitude compensation ensure accurate pressure output and optimal treatment effect.
- With Data Management Software, usage profile (pressure setting, total blower hours, total therapy hours, etc.) can be downloaded via USB interface, allowing physician and homecare provider to track usage record
- Built-in alarm clock

Optional Accessories:



Heated Humidifier



Carry Bag



USB Cable



12V/24V DC to AC Inverter



Specifications:

Pressure Range (cmH2O)	4.0 - 20.0(0.5 cmH2O increment)
Pressure accuracy (cmH2O)	±0.5
Ramp Time (min)	0-45 min(5 min increments)
Ramp Starting Pressure (cmH2O)	3.0 – 19.0
Data Storage	1. Total Meter (Total Blower Hours): 9999.9 hrs 2. Compliance Meter (Total Therapy Hours): 9999.9 hrs
Sound Level	Under 28 dBA*
Dimension (W x D x H)	5.8" x 5" x 3.9"(14.8 cm x 12.8 cm x 10 cm)
Weight	1.76 lb (800 g)
Power Supply	100V-240V, 50/60Hz, Universal power supply

Specifications are subject to change for improvement without notice.

APEX MEDICAL CORP.

9, Min Sheng St. Tu-Cheng,
 Taipei County, 236, Taiwan, R.O.C.
 TEL: 886-2-22683100
 FAX: 886-2-22686525
 E-mail:marketing@apexmedical.com.tw
 www.apexbrand.com

APEX MEDICAL S.L.

Máximo Aguirre 18 Bis, 8ª planta
 48011 Bilbao (Vizcaya)
 Tel: (+34) 94.470.64.08
 Fax: (+34) 94.470.64.09
 E-mail: info@apexmedical.es
 www.apexbrand.com

