

Sino

R104

Happy life
IDEAL

Shopping tray on the seat
Shopping basket in the front of rollator
Powder coating steel frame
Solid 8" x 2" wheels

Seat height:	61 cm
Height adjustable:	79-99 cm
Between handles:	43 cm
Total width:	58 cm
Total length:	71 cm
N.W.:	12.5 kg
Max. weight capacity:	120 kg
Qty./Carton:	1 pc
Carton size:	60 x 27 x 68 cm
Qty./Container:	610 pcs/40'hc



www.idealanycare.com



Pull to fold frame
Release to fix frame



Cane holder



New seat plate