

BABYCROSS



Carrozzina elettrica basculante a telaio fisso da interni ed esterni. Indicata per bambini e ragazzi in età evolutiva in grado di gestire elementi elettronici di comando. Possibilità di personalizzazione attraverso accessori che migliorano il comfort e la sicurezza. Le diverse regolazioni la rendono, inoltre, altamente adattabile al corpo del paziente e ideale per accogliere sistemi posturali complessi. Indicata per terreni particolarmente sconnessi.

Tilting power wheelchair with the fixed frame by internal and external. Suitable for children and teenagers in children can handle electronic control elements. Possibility of customization through accessories that improve comfort and security. The various adjustments make it also highly adaptable to the patient's body and ideal to accommodate complex postural systems. Especially indicated for rough terrain.

Classificazione ISO / ISO classification

12.21.27.009
12.24.03.109
12.24.06.160
12.24.06.163
12.24.06.169
18.09.39.103
18.09.39.133
18.09.39.136

RDM 859599

Utilizzo Use	Pensata per Thought for
Autonomo Independent	Young Young
Outdoor Outdoor	



Caratteristiche base

- basculamento elettronico
- schienale regolabile meccanicamente ed abbattibile per il trasporto
- schienale tensionabile e regolabile in altezza
- seduta regolabile in profondità
- braccioli e fianchetti estraibili e regolabili in altezza ed in larghezza
- pedane estraibili regolabili in altezza e profondità
- indicatori di direzione led

Prestazioni

- velocità max 10km/h
- pendenza massima superabile 15%
- rampe brevi 20%
- autonomia 5h a marcia continua

Motorizzazione

- motori da 300W
- freno elettromagnetico
- regolatore di velocità
- 2 batterie da 12V - 80A

Ruote

- anteriori Ø 280x45 - 8"x2"
- posteriori Ø 350.8 pneumatiche

Base features

- power tilting system
- mechanically adjustable backrest, collapsible for transport
- backrest with tensioning belts, height adjustable
- seat in depth adjustable
- removable armrests and armrest borders adjustable in height and width
- removable footrests, adjustable in height and depth
- led direction indicators

Performance

- max speed 10km/h
- max slope 15%
- short ramps 20%
- autonomy 5h continue march

Motor

- motors 300W
- electromagnetic brake
- speed regulator
- 2 batteries of 12V - 80A

Wheels

- front Ø 280x45 - 8"x2"
- rear Ø 350.8 pneumatic

Accessori & optional / Accessories & optional



Colori disponibili / Available colours



Caratteristiche tecniche / Technical features

