

FROG



Carrozzina elettrica basculante a telaio fisso per gli spazi interni. Indicata per bambini e ragazzi in età evolutiva in grado di gestire elementi elettronici relativi alle azioni di guida. Il sistema di motorizzazione di cui è dotata permette di muoversi in autonomia anche in spazi ridotti. Dotata di sistemi di regolazione multipla per rispondere a esigenze di comfort e affidabilità. La struttura del telaio, inoltre, permette di accogliere sistemi posturali complessi.

Tilting power wheelchair fixed frame to the interior spaces. Suitable for children and in young people can handle electronic elements for driving actions. The operator system which is equipped allows to move autonomously even in small spaces. Equipped with multi-setting systems to meet the requirements of comfort and reliability. The frame structure also allows to accommodate complex postural systems.

Classificazione ISO / ISO classification

12.21.27.006
12.24.03.109
12.24.06.160
12.24.06.163
12.24.06.169
18.09.39.103
18.09.39.133
18.09.39.136

RDM 859598

Utilizzo Use	Pensata per Thought for
Autonomo Independent	Young Young
Indoor Indoor	



Caratteristiche base

- basculamento elettrico
- schienale regolabile meccanicamente ed abbattibile per il trasporto
- schienale tensionabile e regolabile in altezza
- seduta regolabile in profondità
- braccioli e fianchetti estraibili e regolabili in altezza ed in larghezza
- pedana unica regolabile in altezza e flessione
- appoggia gambe imbottito
- cinturini ferma tallone

Prestazioni

- velocità max 8km/h
- pendenza massima superabile 10%
- rampe brevi 15%
- autonomia 3h a marcia continua

Motorizzazione

- motori da 175W
- freno elettromagnetico
- regolatore di velocità
- 2 batterie da 12V - 42"

Ruote

- anteriori Ø 200x45 - 8"x2"
- posteriori Ø 300.7 pneumatiche - 12 1/2"x 2 1/4"

Base features

- power tilting system
- mechanically adjustable backrest and collapsible for transport
- backrest with tensioning belts, in height adjustable
- seat in depth adjustable
- removable armrests and armrest borders, in height and width adjustable
- single footrest adjustable in height and flexion-extension
- padded calf pad
- heel stopping belts

Performance

- max speed 8km/h
- max slope 10%
- short ramps 15%
- autonomy 3h continue march

Motor

- motors of 175W
- electromagnetic brake
- speed regulator
- 2 batteries of 12V - 42A

Wheels

- front Ø 200x45 -8"x2"
- rear Ø 300.7 pneumatic - 12 1/2"x 2 1/4"

Caratteristiche tecniche / Technical features



Accessori & optional / Accessories & optional



Colori disponibili / Available colours

