

PIEGHEVOLI / FOLDABLE

GR108



Carrozzina pieghevole leggera con struttura a doppia crociera. Bracciolidesk estraibili. Pedane girevoli, estraibili e regolabili in altezza. Posizionamento variabile in più posizioni dei perni delle ruote grandi. Posizionamento variabile in più posizioni delle forcelle delle ruote piccole con inclinazione regolabile. Ruota posteriore con perno a destrazione rapida. Ruotine da transito per passaggi stretti.

Lightweight folding wheelchair. Double crossbar structure. Detachable desk size armrests. Swinging, detachable and height-adjustable footrests. Rear wheelshubs adjustable to different positions. Small wheels castors adjustable to different positions with adjustable inclination. Quick-release back wheels. Small wheels for narrow passages.

Classificazione ISO / ISO classification

12.21.06.039
12.24.21.106
12.24.21.109
12.24.21.112
12.24.03.118
12.24.03.121
12.24.06.172

Utilizzo Use	Pensata per Thought for
Autonomo Independent	Senior Senior
Indoor Indoor	
Outdoor Outdoor	



Caratteristiche base

- braccioli estraibili
- pedane estraibili, girevoli e regolabili in altezza
- ruote anteriori piroettanti
- fascia parapalpacci
- freni di stazionamento

Base features

- removable armrests
- removable footrests, swivel and adjustable in height
- twisting front wheels
- calf panel
- service brake

Accessori & optional / Accessories & optional



Spondina desk
Desk armrest



Spondina con bracciolo lungo
Long armrest



Tavolo bordo e incavo
Anatomic tray table



Tavolo rettangolare
Rectangular tray table



Pedana unica
Single footplate



Coppia di pedane elevabili
Pair of swinging
detachable elevating



Sedile con foro WC
Commode seat with inset
cushion



Sedile rigido
Removable padded seat



Spinta tronco destra e
spinta tronco sinistra
Lateral support pad



Cintura semplice
Belt



Cintura di sicurezza
Safety belt



Cintura semplice
con cuscinetto
Padded belt



Bracciolo emiplegico
Armrest for hemiplegic



Appoggiatesta polifunzionale
imbottito
Polyfunctional padded
headrest



Prolunga schienale
Backrest extension



Schienale tensionabile
Tensioning backrest



Sistema monoguida
One-arm driving system



Fascione poggiatesta
Padded leg band



Asta porta flebo
IV pole

Caratteristiche tecniche / Technical features



cm 36 - 38
40 - 43
45 - 48



cm 54 - 56
58 - 61
63 - 66




cm 76



cm 102



cm 41




cm 51



cm 22...30



cm 41



kg 14,5



cm 26
81
92



Optional
cm 45