

GR110



Carrozzina in telaio leggero pieghevole, indicata per bambini e ragazzi in grado di realizzare in modo autonomo o semi autonomo percorsi in ambienti interni ed esterni. Lo schienale regolabile in tensione facilita il contenimento del tronco favorendo maggior comfort. Possibilità di utilizzo in molteplici situazioni grazie alle ruote da transito ed ai perni ad estrazione rapida. Allestita con ruote antiribaltamento che ne migliorano la stabilità.

Folding wheelchair lightweight frame, suitable for children and young people can realize autonomously or semi-autonomous paths in indoor and outdoor environments. The tensioning backrest facilitates the containment of the trunk favouring greater comfort. It can be used in multiple situations thanks to the transit wheels and the quick-release axles. Fitted with antitipping wheels, which improve stability.

Classificazione ISO / ISO classification

12.21.06.036
12.24.03.121
12.24.06.160
12.24.06.172
12.24.21.106
12.24.21.109
12.24.21.112

RDM 858616

Utilizzo Use	Pensata per Thought for
Autonomo Independent	Young Young
Indoor Indoor	
Outdoor Outdoor	



Caratteristiche base

- schienale tensionabile
- regolazione posizione perni delle ruote posteriori
- perni ad estrazione rapida
- forcelle anteriori regolabili in altezza ed inclinazione
- braccioli estraibili e regolabili in altezza
- maniglie di spinta regolabili in altezza
- pedana unica regolabile in altezza, profondità e flesso estensione
- routine antiribaltamento

Base features

- backrest with tensioning belts
- adjustable position of the rear wheels axle
- quick-release axles
- front fork adjustable in height and declivity
- removable armrests, adjustable in height
- pushing handles in height adjustable
- single footplate adjustable in height, depth and flexion-extension
- anti-tipping wheels

Accessori & optional / Accessories & optional



Copriraggi "Tiger" - 22"
Cover rays "Tiger" - 22"



Copriraggi "Spider" - 22"
Cover rays "Spider" - 22"



Appoggiatesta polifunzionale
Polyfunctional headrest



Appoggiatesta polifunzionale in poliuretano piccolo
Polyfunctional PU small size headrest



Cintura semplice
Belt



Cintura semplice con cuscinetto
Padded belt



Cintura pelvica piccola
Small size pelvic belt



Spinta tronco destra e sinistra
Lateral support pad



Tavolo in plexyglass
Plexyglass tray table

Colori disponibili / Available colours



RAL 9001 RAL 1028 RAL 2004 RAL 3020 RAL 4005 RAL 4010 RAL 5002 RAL 5012 RAL 6018

Caratteristiche tecniche / Technical features



cm 23 33 cm 33 cm 36 cm 49 cm 48 56 cm 80 cm 17 - 25 kg 15 kg 60
28 36
30 59
53