

GR120



Carrozzina superleggera e affidabile, con struttura pieghevole in alluminio verniciato a doppia crociera, disponibile in vari colori. Dotata di schienale tensionabile, forcelle regolabili e perni ad estrazione rapida. Presenta caratteristiche tecniche che permettono la personalizzazione dell'assetto ergonomico al fine di ottimizzarne comfort e stabilità. La sua struttura flessibile è facile da manovrare grazie al peso limitato e alle ruote di transito di cui è dotata.

Superlight and reliable wheelchair, folding structure in painted aluminium, with double cross, available in various colors. Tensioning backrest equipped with adjustable forks and quick-release axles. It presents technical features that enable ergonomic trim customization in order to maximize comfort and stability. Its flexible structure is easy to manoeuvre thanks to low weight and transit wheels supplied with it.

Classificazione ISO / ISO classification

12.21.06.060

RDM 858621

Utilizzo Use	Pensata per Thought for
Autonomo Independent	Senior Senior
Indoor Indoor	
Outdoor Outdoor	




Caratteristiche base

- telaio in alluminio (alu 5000)
- doppia crociera
- schienale tensionabile
- regolazione posizione perni delle ruote posteriori
- perni ad estrazione rapida
- forcelle anteriori regolabili in altezza ed inclinazione
- spondine para abiti
- pedane fisse regolabili in altezza
- ruote da transito


Base features

- aluminum structure (alu 5000)
- double crossbar
- backrest with tensioning belts
- adjustable position of the rear wheels axle
- quick-release axles
- front forks adjustable in heights and declivity
- sideboards protecting clothes
- fixed footrests adjustable in height
- small wheels for narrow passages

Accessori & optional / Accessories & optional



Colori disponibili / Available colours



Caratteristiche tecniche / Technical features

