

DOCCIA / SHOWER

GR312



Colori disponibili / Available colours



RAL 9001

RAL 6027

Accessori & optional / Accessories & optional

Cintura semplice
BeltCintura pelvica grande
Big size pelvis belt

Caratteristiche tecniche / Technical features



cm 45



cm 43



cm 48



cm 60



cm 58



cm 86



cm 24



kg 22



kg 120

Sedia da doccia basculante per persone affette da grave disabilità motoria che necessitano di assistenza durante le azioni di trasferimento ed igiene personale.

Tilting shower chair for people affected by serious movement disability, who needs assistance during transfer and personal hygiene.

Classificazione ISO / ISO classification

09.12.03.003

18.09.18.202

18.09.18.205

RDM 858734

Utilizzo
UseFascia d'età
Age groupAssistito
AssistedSenior
SeniorIndoor
IndoorSilver
Silver

Caratteristiche base

- sistema basculante
- ruote piroettanti
- dispositivi di stazionamento su tutte le ruote
- sedile prodotto con materiale idrorepellente ed ammortizzante
- schienale in tela tensionabile
- appoggiatesta
- braccioli ribaltabili
- pedane regolabili in altezza
- fascia appoggia gambe
- vaschetta WC estraibile

Base features

- tilting system
- twist wheels
- service brake on every wheels
- seatrest produced with waterproof and cushioning material
- canvas backrest with belts
- headrest
- turning up armrests
- footrests, adjustable in height
- legs leaning belt
- removable WC bowl

MADE IN
ITALY