

GR205






Poltrona imbottita dotata di due ruote piroettanti per uso assistito. Indicata per persone che necessitano dell'assistenza di un accompagnatore in ambienti interni ed esterni. Ideale per spostamenti in spazi ridotti. Dotata di schienale avvolgente, sagomato anatomicamente, per accogliere il tronco e prevenire situazioni di instabilità. Possibilità di personalizzazione attraverso una vasta gamma di accessori che ne migliorano funzionalità e sicurezza.

Padded chair with two swivel wheels for assisted use. Indicated for people who need the assistance of an attendant for indoor and outdoor environments. Ideal for moving around in tight spaces. It is equipped with comfortable backrest anatomically shaped to accommodate the trunk and prevent situations of instability. Possibility of customization through a wide range of accessories that improve functionality and security.

Classificazione ISO / ISO classification

12.21.03.006
12.24.21.103

RDM 858686

Utilizzo Use	Pensata per Thought for
Assistito Assisted 	Senior Senior
Indoor Indoor 	Silver Silver
Outdoor Outdoor 	

MADE IN
ITALY 

Caratteristiche base

- schienale imbottito avvolgente e reclinabile
- braccioli estraibili
- sedile imbottito con foro per WC
- vaschetta WC estraibile
- pedane elevabili e regolabili in altezza
- poggiatesta imbottiti
- ruote anteriori Ø 200mm
- ruote posteriori piroettanti Ø 175mm

Base features

- padded enveloping and reclining backrest
- removable armrests
- padded seat with WC hole
- removable WC bowl
- elevating footrests and adjustable in height
- padded legrests
- front wheels Ø 200mm
- rear twisting wheels Ø 175mm

Accessori & optional / Accessories & optional



Appoggiatesta imbottito
Padded headrest

Bracciolo emiplegico
Armrest for hemiplegic

Pedana unica
Single footplate

Sedile con foro a "U"
Seat with a "U" hole

Cintura semplice
Belt

Cintura pelvica grande
Big size pelvis belt

Spinta tronco destra e
spinta tronco sinistra
Lateral support pad

Cuneo divaricatore
imbottito
Padded abductor pad

Tavolo bordo e incavo
Anatomic tray table

Tavolo imbottito
Padded table

Tavolo rettangolare
Rectangular tray table


Asta portaflebo
IV drip stand

Colori disponibili / Available colours



RAL 7015

Caratteristiche tecniche / Technical features



cm 43
48
53

cm 44

cm 90

kg 120