



TOP VIEW



SIDE VIEW



Un-Controlled Copy

.TD.

1. BLANK FORCEPS : PART # 20-11000-83 : STAINLESS STEEL : SS410/420
2. PLASTIC CAP : PC 105 : MOLDED PLASTIC BLACK : DELRIN OR EQUIVILANT
3. CAP FILLING : HI-TEMP EPOXY : BLACK : EPON RESIN 828 OR EQUIVILANT
4. INSULATION : NYLON : (PVDF) BLUE OR AS PER SPECIAL INSTRUCTIONS

TOLERANCES ON		INCHES
DECIMALS	XXX +/- .005	
	XX +/- .015	+/- .25
	X +/- .030	
OR AS STATED ON DWG.		
DRAWN BY: <i>[Signature]</i>	Date	4/15/97
CHECKED BY: <i>[Signature]</i>	Date	4/16/97
ENGINEER (APPROVED) <i>[Signature]</i>		



ARTICOLI & APPARECCHI PER MEDICINA

6" STRAIGHT FORCEPS

REV #	ECN #	DRAWING No.	SCALE
N/A	N/A	310-120	1" = 1"