

Pieghevole per facile trasporto in auto

Sedie elettroniche Electronic wheelchairs



Box batterie con rotelle per facile trasporto
Batteries, located in a box with wheels for easy moving



700 EXTRA

Sedia elettronica pieghevole per esterno

Standard: Motori da 330 Watt con batterie da 100 A senza manutenzione. Comando con clacson, regolatore di velocità. Telaio pieghevole. Leggere imbottiture con tessuti ignifughi. Freni elettromagnetici. Braccioli regolabili in altezza ed estraibili. Fanaleria regolamentare con indicatori di direzione. Cintura di sicurezza. Ruote antiribaltamento. Caricabatteria. Velocità max 8 km/h. Autonomia 6h - 40 km.

Optional: Appoggiambraccioli regolabili, girevoli ed estraibili. Schienale abbattibile. Batterie al litio, +60% di durata, -18 Kg di peso.

700 EXTRA

Electronic folding chair for external use

Standard: 330 Watt motors with 100 A batteries without maintenance. Control box with horn, speed governor. Folding frame. Light padding with fire-retardant upholsteries. Electromagnetic brakes. Adjustable in height and removable armrests. Regulatory lights and indicators. Security belt. Anti-overflow small wheels. Battery charger. Max speed 8 km/h. Range: 6h - 40 km.

Optional: Elevating, swinging and detachable legrests. Foldable backrest. Lithium batteries, life 60% more, weight 18 Kg less.

Codice LEA 12.36.06.015

ISO Class 12.21.27.009 ricond.+12.24.03.163

