

Nuova **ITALA** by **SURACE**

La sedia leggera di qualità
ad un prezzo imbattibile

Caratteristiche Itala

- Telaio in alluminio ed acciaio speciale ad alta resistenza
- Monocrociera rinforzata
- Braccioli desk girevoli
- Pedane regolabili, girevoli ed estraibili
- Regolazione del sedile in altezza

Optional

- Kit speed (passo accorciato e rotelline antiribaltamento)
- Ruote ad estrazione rapida con ruote ascensore (di serie su mod. SUPERITALA)
- Catarifrangenti
- Cintura di sicurezza




Ruote ad estrazione rapida

...e **SUPERITALA** con ruote pneumatiche ad estrazione rapida e ruote ascensore, in versione antracite di serie.

15 kg



ITALA Class. ISO: 12.21.06.039 + 12.24.03.121 SUPERITALA Class. ISO: 12.21.06.039 + 12.24.03.121 + 12.24.21.106



SURACE S.p.A.
C.so Savona 52, 10024 Moncalieri (TO) - Italy
Tel: +39 011 6408476 - Fax: +39 011 641562
surace@surace.it - www.surace.it

SURACE
WHEELCHAIRS

New **ITALA** by **SURACE**

The quality lightweight wheelchair
at an unbeatable price.

Itala features:

- aluminium and special high-strength steel frame
- single reinforced cross
- swinging and desk armrests
- swinging, adjustable and removable footrests
- adjustable in height seat

Optional:

- speed kit (shortened step with anti-overturn small wheels)
- extractable big wheels and wheels for lift (standard on SUPERITALA)
- reflectors
- security belt


... and **SUPERITALA** with pneumatic, extractable big wheels and wheels for lift, with standard frame in dark grey colour

15 kg



Extractable big wheels

Itala ISO class: 12.21.06.039 + 12.24.03.121 Superitala ISO class: 12.21.06.039 + 12.24.03.121 + 12.24.21.106



SURACE S.p.A.
C.so Savona 52, 10024 Moncalieri (TO) - Italy
Tel: +39 011 6408476 - Fax: +39 011 641562
surace@surace.it - www.surace.it

SURACE
WHEELCHAIRS