

Classica e bella per ogni uso  
Larghezza esterna solo 58 cm



Sedie elettroniche Electronic wheelchairs



=



+ Kg 30  
BATTERY

**710 MAGIC**  
Sedia elettronica pieghevole universale

**Standard:** Motori da 330 Watt con batterie da 100 A senza manutenzione. Comando con clacson, regolatore di velocità. Telaio pieghevole. Fianchi desk girevoli. Appoggiamambe regolabili, girevoli ed estraibili. Leggere imbottiture con tessuti ignifughi. Freni elettromagnetici. Cintura di sicurezza. Ruote antiribaltamento. Caricabatteria. Velocità max 7 km/h. Autonomia 6h - 40 km.

**Optional:** Comando retrattile. Comando a mento. Comando allo schienale. Schienale abbattibile.

**710 MAGIC**  
Electronic folding universal chair

**Standard:** 330 Watt motors with 100 A batteries without maintenance. Control box with horn, speed governor. Folding frame. Desk and tipping armrests. Adjustable, swinging and detachable legrests. Light padding with fire-retardant upholsteries. Electromagnetic brakes. Security belt. Anti-overturn small wheels. Battery charger. Max speed 7 km/h. Range: 6h - 40 km.

**Optional:** Retractable control box. Chin control box. Control box on backrest. Foldable backrest.

Codice LEA 12.23.06.009

ISO Class 12.21.27.006

cm Ø 33	cm Ø 21	L 44	L 58	L 75 - 108	h 52	Standard	20%	Optional	Optional	Optional	Optional	Optional	Portata