

**Nome Prodotto**

Barella per Sollevatore Duo.

Codice	Misure (cm)	Portata (Kg)
01.01.1100	185 x 55 x 10	160

**Destinazione d'uso**

Le Barelle vengono associate ad un sollevatore della Linea SCANMEDICALEUROPE per la movimentazione dei pazienti.

**Descrizione**

La BARELLA DUO a cucchiaio ha un telaio in acciaio con doppia chiusura a baionetta e lastre d'appoggio in plexiglass. La barella è radiotrasparente.

Materiali	
Telaio	Acciaio
Lastre d' Appoggio	Plexiglass
Manutenzione e Pulizia	
Disinfezione possibile con disinfettanti diluiti al 5% anche a base alcolica	
Conservazione e Smaltimento	
Effettuare gli spostamenti solamente dopo aver controllato che il paziente sia stato agganciato in modo corretto e sicuro (vedi manuale d'istruzioni)	
Conservare al riparo da fiamme libere e fonti di calore	
Smaltire il dispositivo rispettando le normative vigenti in materia di tutela dell'ambiente e raccolta differenziata	
Confezione e Nr. Articoli	
Scatola	1 Pezzo
Garanzia	
Europea	24 mesi per difetti di fabbricazione ad esclusione delle parti consumabili
Riferimenti del Prodotto	
Fabbricante	Demarta Virginio S.a.S.
Classe di Dispositivo 93/42	Classe I
Codice CND	Y123603
Codice RDM	1689778/R

### Product Name

Stretcher Duo.

Code	Dimensions (cm)	Capacity (Kg)
01.01.1100	185 x 55 x 10	160



### Use

The stretcher are associated with a lift of the SCANMEDICALEUROPE Line for patient handling.

### Description

The spoon STRETCHER DUO has a steel frame with double bayonet lock and plexiglass support plates.  
The stretcher is radiolucent.

Materials	
Frame	Steel
Support Plates	Plexiglass
Maintenance and Cleaning	
Disinfection is possible with 5% diluted disinfectants even with alcohol base	
Storage and Disposal	
Make the movements only after checking that the patient has been correctly and securely locked (see instruction manual)	
Keep away from heat, sparks, open flame or other sources of ignition	
Dispose of the device in compliance with current regulations on environmental protection and separate waste collection	
Package and Nr.Items	
Box	1 Piece
Warranty	
European	24 months for manufacturing defects with the exception of consumable parts
Product References	
Manufacturer	Demarta Virginio S.a.S.
Directive 93/42	Class I
CND Code	Y123603
RDM Code	1689778/R