

Nome Prodotto

Carrozzina ad Autospinta Viky Exchange Plus.

Destinazione d'uso

Le carrozzine sono un mezzo di trasporto indispensabile per le persone disabili con assente o ridotta capacità di deambulazione.

Lo scopo è quello di facilitare gli spostamenti autonomi e di salvaguardare il benessere fisico della persona, evitando la comparsa di complicanze correlate alla posizione seduta obbligatoria (piaghe da decubito, dolori posturali cronici...).



Codice	Taglia	Lunghezza	Larghezza	Altezza	Larghezza Piegata	Larghezza Seduta	Profondità Seduta	Altezza Schienale	Altezza Seduta da Terra	Portata Massima (Kg)	Dimensioni Ruote Anteriori	Dimensioni Ruote Posteriori	Peso Netto con Poggiatesta	Dimensioni Scatola
Misure in Cm														
01.04.41EDQR41FD	41	80	56	90.5	26	41	40	44	50	120	20 (8")	60 (24")	18.7	82x27x93
01.04.41EDQR43FD	43	80	58.5	90.5	26.5	43	40	44	50	120	20 (8")	60 (24")	18.8	82x27x93
01.04.41EDQR46FD	46	80	61	90.5	27	46	40	44	50	120	20 (8")	60 (24")	18.9	82x27x93
01.04.41EDQR48FD	48	80	64	90.5	27.5	48	40	44	50	120	20 (8")	60 (24")	19	82x27x93
01.04.41EDQR50FD	50	80	65	90.5	28	50	40	44	50	120	20 (8")	60 (24")	20	82x27x93

Caratteristiche

Telaio	Acciaio Verniciato a Polvere
Crociera	Doppia
Colore	Grigio Canna di Fucile (Opaco)
Caratteristiche del Telaio	Doppio foro sul telaio per la trasformazione della carrozzina da autospinta in transito
Ruote Anteriori	20 x 2,5 cm (8" x 1") PVC Grigie
Ruote Posteriori	Estraibili da 60 cm (24") PU Grigie
Accessori	Ruote Anti-ribaltamento e Ruotine per Ascensore
Braccioli	Modello a Scrivania, Ribaltabili
Pedane	Estraibili e Regolabili con Supporto per Gambe
Poggiatesta	PVC Nero
Sedile e Schienale	Nylon Morbido Imbottito, lavabile e disinfettabile

Manutenzione e Pulizia

Disinfezione possibile con disinfettanti diluiti al 5% anche a base alcolica

Conservazione e Smaltimento

Si consiglia l'uso ad un solo utente

Conservare al riparo da fiamme libere e fonti di calore

Smaltire il dispositivo rispettando le normative vigenti in materia di tutela dell'ambiente e raccolta differenziata

Confezione e Nr. Articoli

Scatola	1 Pezzo
---------	---------

Garanzia

Europea	24 mesi per difetti di fabbricazione ad esclusione delle parti consumabili
---------	--

Riferimenti del Prodotto

Fabbricante	Demarta Virginio S.a.S.
Classe di Dispositivo 93/42	Classe I
Codice ISO	12.22.03.006
Codice CND	Y122106
Codice RDM	1721026/R

Product Name

Manual Wheelchairs Viky Exchange Plus.

Use

The wheelchairs are an essential way of transport for disabled people with none or limited walking ability.

The purpose is to facilitate autonomous movements and to safeguard the physical well-being of the person, avoiding the appearance of complications related to the sitting position (bed sores, chronic postural pains ...).



Code	Size	Length	Width	Height	Width Folded	Seat Width	Seat Depth	Back Height	Floor to Seat Height	Weight Capacity	Front Wheel	Rear Wheel	Weight with Footrest	Carton Size
Measures in Cm														
01.04.41EDQR41FD	41	80	56	90.5	26	41	40	44	50	120	20 (8")	60 (24")	18.7	82x27x93
01.04.41EDQR43FD	43	80	58.5	90.5	26.5	43	40	44	50	120	20 (8")	60 (24")	18.8	82x27x93
01.04.41EDQR46FD	46	80	61	90.5	27	46	40	44	50	120	20 (8")	60 (24")	18.9	82x27x93
01.04.41EDQR48FD	48	80	64	90.5	27.5	48	40	44	50	120	20 (8")	60 (24")	19	82x27x93
01.04.41EDQR50FD	50	80	65	90.5	28	50	40	44	50	120	20 (8")	60 (24")	20	82x27x93

Features

Frame	Powder Coating Steel Frame
Cross Bar	Double Cross Brace
Color	Smokey Grey (Dull Polish)
Frame Features	Double hole on the frame for the transformation from manual wheelchair to a transit one
Front Wheels	20 x 2,5 cm (8" x 1") PVC Grey Tyre
Rear Wheels	60 cm (24"), Spoke Wheels, PU Grey Tyre, Quick Release
Accessory	Small Wheels for Lift and Anti Tipper Wheels
Armrest	Flip Back Slope Model Armrest
Footrest	Swing In and Away, Length Adjustable with Leg Support
Footplate	Plastic Black
Seat and Backrest	Soft and Padded Nylon, washable and disinfectable

Maintenance and Cleaning

Disinfection is possible with 5% diluted disinfectants even with alcohol base

Storage and Disposal

Single-user Product

Keep away from heat, sparks, open flame or other sources of ignition

Dispose of the device in compliance with current regulations on environmental protection and separate waste collection

Package and Nr. Items

Box 1 Piece

Warranty

European 24 months for manufacturing defects with the exception of consumable parts

Product References

Manufacturer	Demarta Virginio S.a.S.
Directive 93/42	Class I
ISO Code	12.22.03.006
CND Code	Y122106
RDM Code	1721026/R