

Nome Prodotto















Carrozzina ad Autospinta Viky Slim.

Destinazione d'uso

Le carrozzine sono un mezzo di trasporto indispensabile per le persone disabili con assente o ridotta capacità di deambulazione.

Lo scopo è quello di facilitare gli spostamenti autonomi e di salvaguardare il benessere fisico della persona, evitando la comparsa di complicanze correlate alla posizione seduta obbligatoria (piaghe da decubito, dolori posturali cronici...).



Codice	Taglia	Lunghezza	Larghezza	Altezza	Larghezza Piegata	Larghezza Seduta	Profondità Seduta	Altezza Schienale	Altezza Seduta da Terra	Portata Massima (Kg)	Dimensioni Ruote Anteriori	Dimensioni Ruote Posteriori	Peso Netto con Poggiatesta	Dimensioni Scatola
Misure in Cm														
01.04.340D41FD	41	69	51	97	30.5	41	40	40	54	120	20 (8")	51 (20")	18.5	75x34x99.5
01.04.340D43FD	43	69	53	97	31	43	40	40	54	120	20 (8")	51 (20")	18.6	75x34x99.5
01.04.340D46FD	46	69	55	97	31.5	46	40	40	54	120	20 (8")	51 (20")	18.7	75x34x99.5

Caratteristiche

Telaio	Acciaio Verniciato a Polvere
Crociera	Doppia
Colore	Blu-Nero
Ruote Anteriori	20 x 2,5 cm (8" x 1") PVC Grigie
Ruote Posteriori	51 cm (20") PU Grigie
Anello Corrimano	Acciaio Verniciato a Polvere Nero
Braccioli	Estraibili, Modello a Scrivania con Poggiabraccio in PVC
Pedane	Estraibili e Regolabili con Supporto per Gambe
Poggiatesta	PVC Nero
Sedile e Schienale	Nylon Morbido Imbottito, lavabile e disinfettabile

Manutenzione e Pulizia

Disinfezione possibile con disinfettanti diluiti al 5% anche a base alcolica

Conservazione e Smaltimento

Si consiglia l'uso ad un solo utente

Conservare al riparo da fiamme libere e fonti di calore

Smaltire il dispositivo rispettando le normative vigenti in materia di tutela dell'ambiente e raccolta differenziata

Confezione e Nr. Articoli

Scatola 1 Pezzo

Garanzia

Europea 24 mesi per difetti di fabbricazione ad esclusione delle parti consumabili

Riferimenti del Prodotto

Fabbricante	Demarta Virginio S.a.S.
Classe di Dispositivo 93/42	Classe I
Codice ISO	12.22.03.003
Codice CND	Y122106

Product Name


















Manual Wheelchairs Viky Slim.

Use

The wheelchairs are an essential way of transport for disabled people with none or limited walking ability.

The purpose is to facilitate autonomous movements and to safeguard the physical well-being of the person, avoiding the appearance of complications related to the sitting position (bed sores, chronic postural pains ...).



Code	Size	Length	Width	Height	Width Folded	Seat Width	Seat Depth	Back Height	Floor to Seat Height	Weight Capacity	Front Wheel	Rear Wheel	Weight with Footrest	Carton Size			
Measures in Cm																	
01.04.340D41FD	41	69	51	97	30.5	41	40	40	54	120	20 (8")	51 (20")	18.5	75x34x99.5			
01.04.340D43FD	43	69	53	97	31	43	40	40	54	120	20 (8")	51 (20")	18.6	75x34x99.5			
01.04.340D46FD	46	69	55	97	31.5	46	40	40	54	120	20 (8")	51 (20")	18.7	75x34x99.5			

Features	
Frame	Powder Coating Steel Frame
Cross Bar	Double Cross Brace
Color	Blue-Black
Front Wheels	20 x 2,5 cm (8" x 1") PVC Grey Tyre
Rear Wheels	51 cm (20"), Spoke Wheels, PU Grey Tyre
Handrim	Powder Coating Black Steel
Armrest	Removable Armrest with PVC Airmat
Footrest	Swing Away and Length Adjustable with Leg Support
Footplate	Plastic Black
Seat and Backrest	Soft and Padded Nylon, washable and disinfectable
Maintenance and Cleaning	
Disinfection is possible with 5% diluted disinfectants even with alcohol base	
Storage and Disposal	
Single-user Product	
Keep away from heat, sparks, open flame or other sources of ignition	
Dispose of the device in compliance with current regulations on environmental protection and separate waste collection	
Package and Nr. Items	
Box	1 Piece
Warranty	
European	24 months for manufacturing defects with the exception of consumable parts
Product References	
Manufacturer	Demarta Virginio S.a.S.
Directive 93/42	Class I
ISO Code	12.22.03.003
CND Code	Y122106