

La migliore soluzione per il trasporto su aeromobile, treno, ambulanza
The best solution for transport on aircraft, train, ambulance

Con sedile
profondo
e confortevole
*With deep and
comfortable
seat*



Air 2000 Surace: solo 35 cm di ingombro massimo.
Air 2000 Surace: only 35 cm of width.



Barre portantina (optional)
Sedan bars (optional)

AIR 2000
Sedia per il transito sull'aeromobile

Standard: Telaio in acciaio rinforzato. Prolunga appoggiatesta facilmente piegabile. Ruote anteriori piroettanti. 2 cinture di sicurezza allo schienale. 1 cintura di sicurezza al sedile. Fascia ai talloni. Tessuti ignifughi.

Optional: Freni a tamburo. Cassettino sotto sedile porta documenti. Barre a portantina. Logo cliente retro schienale.

AIR 2000
Chair for transit on the aircraft

Standard: Reinforced steel frame. Easily foldable headrest extension. Front swivel wheels. 2 security belts on the back. 1 security belt on the seat. Heels strap. Fire proof upholstery.

Optional: Drum brakes. Storage drawer under the seat for documents. Sedan bars. Customer logo at the back.

cm Ø 21	cm Ø 8	h 60 - 90	L 35	h 52	L 40	L 60	Optional	Weight 12 Kg	Portata 130 kg